



TCP Wizard 96+ Hour Frozen or Refrigerated Shipping Solution




- Frozen or Refrigerated/-16° or 0°C
- For up to 96 hours
- No dry ice or wet ice
- Designed for multiple use
- Used by clinical research studies, medical & biomedical laboratories

IATA approval pending

Payload Space: 9.5" x 9.5" x 9.5"

The TCP Wizard, maintains desired payload temperature conditions for shipments of laboratory specimens over long distances which require prolonged transport times of one to four days. Internal temperature keeps the laboratory specimens integrity regardless of external ambient environment. The Wizard requires no dry ice or messy wet ice. Additionally, our proprietary, USA made phase change materials can be easily frozen within 24 hours inside any standard freezer. The Wizard is reusable and contains no hazardous materials thereby eliminating any hazmat training for use and requirements for a hazmat disclosure statement for commercial shipments. For domestic use, refrigerated and frozen shipments can be sent second day air or three day ground, instead of overnight at a considerable less cost for shipping.



The **black** graph line represents the recorded outside temperatures during transit.

The **red** graph line represents the recorded inside payload temperature during transit.